

Breathe-Link K5 Change Logo Procedure

1. Overview

Breathe-Link K5 1.1 software allows users to add their logos on all PDF and print reports. This document outlines procedures to add a logo onto the PDF or print reports





3cm

12.5cm

BreatheLink Session Report

Personal Details:

ID:	R-huygyX-210612-155943	Age:	14
Last Name:		Weight:	145.00
First Name:	huygy	Height:	71.00
Date of Birth:	01/03/98	BMI:	20.22
Gender:	Male	PNV (cmH2O):	93.00

Session Summary:

Date:	21/06/12
Time:	15:59:43
Breaths Completed:	30
Target Load (cmH2O):	67
Session Type:	Auto / Light

2. Procedures

2.1 Logo graphic format

The logo should adhere to the following specification for the best results:

Format type: jpeg
 Size: 3cm height by 12.5cm length
 Name it: your logo

Please note that the software will, if required resize the graphic to fit the space available on the printout.

2.2 Windows XP or later

Locate the BreatheLink folder on the hard drive, normally C:\ drive:

C:\Program Files\BreatheLink\resources

In the **resources** folder replace the existing ' **your logo.jpeg** ' file with the logo you created above. If you are dragging your own logo directly and it is correctly named as ' **your logo.jpeg** ' windows will ask you whether you want to replace existing file. Click yes. Congratulations you have succeeded!

To test whether you have been successful open the Breathe-Link software and generate a report and ensure that the logo appears.

2.3 Windows Vista and Windows 7

Vista and Windows 7 file structure is slightly different to Windows XP or later. You will find the **resources** folder in:

C:\Program Files (x86)\BreatheLink\resources

In the **resources** folder replace the existing ' **your logo.jpeg** ' file with the logo you created above. If you are dragging your own logo directly and it is correctly named as ' **your logo.jpeg** ' windows will ask you whether you want to replace existing file. Click yes. Congratulations you have succeeded!

To test whether you have been successful open the Breathe-Link software and generate a report and ensure that the logo appears.

HINT:

If you are struggling to locate the Breathe-Link folder, locate your Breathe-Link icon on your desktop highlight and right click. Submenu will appear, click on properties. In the location attribute you will see the Breathe-Link path i.e. C:\Program Files (x86)\BreatheLink. Click on the Find Target button, which will navigate you automatically to the BreatheLink folder where the **resources** folder resides. Repeat the above procedure.